



PARTNERS
FOR ENDOSCOPY

Protein Pens



Protein Pens are a simple test with effective results for the detection of protein. An incubator is also available for testing. After an incubation time of at least 45 minutes, the swab in the fluid is examined for the presence of protein or amino acid.

Advantages:

- Portable and convenient
- Shelf stable
- No added chemicals required
- Easy to distinguish positive and negative results

Protein Pens Product Code: SEN3023

Incubator Product Code: SEN3061

For more information please call **01782 594164**
or via email **sales@pfmedical.co.uk**

Hartley House, Galveston Grove, Fenton, Stoke-on-Trent ST4 3PE www.pfmedical.co.uk

All prices are valid for 30 days and exclude VAT and Carriage

Last updated June 2018